

## Information Sheet

### Connection of PR2000S to IA9200, IA9300 Instruments

If you require to use the PR2000S printer on an early model IA9000 instrument (i.e. pre 2012 manufacture), then an alternative connection cable AZ9225 will be needed. This part is available as an accessory and may be ordered via your local distributor for BSL equipment.

Early model IA9000 instruments can easily be identified as they have the model type number ( i.e. 9100, 9200 or 9300 shown on the BLUE KEYPAD. E.g.



(Later models show only the Electrothermal brand name on the keypad. )

For early model IA9000 instruments, refer to Instruction Book M6346 for initial printer set up instructions. If you do not have the M6346 Instruction Book originally supplied with your early model IA, then a PDF copy can be obtained via your local distributor for BSL equipment.